


ThinPrep® UroCyte™

Quick reference guide


(For use with urine cytology processing or slide-based molecular testing)

In order to perform UroCyte sample preparation on the ThinPrep 2000 processor, additional software must be obtained. Please refer to the ordering information section of the manual or call Customer Service for more information.




1. Collection.

Collect urine directly into the ThinPrep UroCyte urine collection kit,* **OR** process urine fresh. (Note: Fresh urine can be mixed with a 2:1 urine to PreservCyt® solution ratio and stored for up to 48 hours before processing.)



5. Evaluate cell pellet appearance.

If cell pellet is not free of blood, add 30ml of CytoLyt solution and repeat from step 4.




2. Concentrate by centrifugation.

- Transfer the sample evenly into two labeled 50ml centrifuge tubes.
- Centrifuge at 600g for 10 minutes or 1200g for 5 minutes.




6. Add the contents of a fresh PreservCyt solution vial (20ml) into the tube.

Vortex briefly to mix the solution and pour the entire sample back into the PreservCyt solution vial. Allow to stand in PreservCyt solution for 15 minutes.




3. Pour off supernatant and resuspend cell pellet.

Resuspension can be done on a vortexor or may be achieved by syringing the pellet back and forth with a plastic pipette.



7. Run on ThinPrep 2000 processor using sequence 5 (UroCyte) or ThinPrep 5000 on sequence UroCyte.

Fix, stain, and evaluate cytology, **OR** perform molecular diagnostic testing according to the manufacturer's instructions for use.




4. Add CytoLyt® solution wash.

- Add 30ml of CytoLyt solution to one 50ml centrifuge tube and vortex. Transfer the contents of this tube into the second 50ml centrifuge tube and vortex. The specimen is now combined into one 50ml tube. The empty tube can be discarded.
- Centrifuge.
- Pour off supernatant.
- Resuspend cell pellet.

*If using the UroCyte collection kit, do not exceed a 2:1 ratio of urine to PreservCyt solution. If the urine volume collected exceeds 60ml, pour off excess. A minimum volume of 33ml of urine is required to perform the UroVysion® assay.


Hologic provides this ThinPrep UroCyte quick reference guide ("guide") as a general informational tool only and is not an affirmative instruction or guarantee. While the information provided in this Guide may describe a particular technique or protocol, it is not intended as an endorsement and/or requirement by Hologic to utilize such technique or protocol. It is the sole responsibility of the laboratory to determine which specific technique and/or protocol to employ. Laboratories using ThinPrep or PreservCyt products must read and understand each product's Instructions for Use and comply with applicable local, state and federal rules and regulations.

Instructions for using the ThinPrep® UroCyte™ urine collection kit



Step 1.

On the specimen collection cup, record patient information in the space provided.



Step 4.

Tightly secure blue cap on specimen cup to prevent leakage. (Keep turning for another 1/4 inch after you hear the audible click.)



Step 2.

Collect urine in a routine manner. If urine volume exceeds 60ml, pour off excess. The total volume of urine must not exceed 60ml. A minimum of 33ml of urine is required to perform the UroVysion® assay.



Step 5.

Place cup and absorbent pads into biohazard bag. Tightly seal bag.



Step 3.

After the urine is collected, carefully pour PreservCyt® solution into specimen cup containing urine. Do not spill PreservCyt solution.

Step 6.

Store between 4°C and 30°C (39°F–86°F). Preferred storage and shipping conditions are on ice packs (e.g., blue ice in Styrofoam™). Specimen must be processed within 48 hours. Transport the specimen according to your internal procedures.

Please refer to ThinPrep, UroCyte, and PreservCyt instructions for complete instructions and information on limitations, warning, and precautions associated with the use of PreservCyt solution.

Diagnostic Solutions | hologic.com | diagnostic.solutions@hologic.com | +1.781.999.7300

For technical support call +1.800.442.9892 - Option 6.

86369-001 Rev. 002 © 2014 Hologic, Inc. All rights reserved. Specifications are subject to change without prior notice. Hologic, CytoLyt, PreservCyt, ThinPrep, UroCyte and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your Hologic representative or write to diagnostic.solutions@hologic.com.